

Q Baby™ Application Instructions

Bio Skin® should be applied to dry skin. Oils, lotions and perspiration will make the support harder to put on and will make it will feel as if it's sized too small. Dry the area of application before applying support. Take the support with both hands and stretch the top and bottom of the support two to three times. This helps to relax the material before its gets broken-in.



1. Position silicone buttress over patellar tendon. Wrap ends, securing the hook on top of opposite end, and secure to a snug fit.



2. Attach the strap to the support with the buckle on either side of the silicone pad. Wrap the elastic strap around the back of the leg. Insert the strap through the buckle and tighten until comfortably snug. Then attach the Velcro®.



3. The elastic strap is tightened to comfort level.